



The Physician Driven Imaging Solution®

Network Connection Guidelines

IT DEPARTMENT MUST ALLOW THE FOLLOWING EXECUTABLE FILES TO RUN:

1. Windows: **Surgimap for Window.exe and Application\Surgimap.exe**
2. MAC: **Surgimap.app**

MUST ALLOW ACCESS TO THE FOLLOWING DOMAINS:

1. <https://www.surgimap.com/api>
2. If the above does not work, please try <http://www.surgimap.com/api>

BELOW ARE THE SPECIFIED SURGIMAP PORTS NEEDED FOR AUTHENTICATION & SYNCHRONIZATION, WEB FILES & UPDATES, AND DOWNLOADING ENCRYPTED IMAGE DATA (MAY BE NECESSARY IF RUNNING SURGIMAP BEHIND HOSPITAL FIREWALL):

1. api.surgimap.com TCP port 8080 SSL
2. www.surgimap.com TCP port 80 HTTP
3. www.surgimap.com TCP port 443 SSL
4. surgimapaccess.s3.amazonaws.com TCP port 80 HTTP
5. surgimapaccess.s3.amazonaws.com TCP port 443 SSL

FOR OUR FDA CLEARANCE, PLEASE SEE BELOW:

<https://goo.gl/zkDIWM>

FOR ANY OTHER QUESTIONS OR INQUIRIES REGARDING THE NETWORK CONNECTION, FIREWALL, AND LOCAL SETUP, PLEASE EMAIL US AT SUPPORT@SURGIMAP.COM.